



Ear Tag Applicator

Weiyou® is a professional ear tag applicator manufacturer and supplier in China. We have been specialized in veterinary products for many years. Weiyou® ear tag applicator have a low-price advantage and cover most of the European and American markets. Weiyou® hope to cooperate with more customers for mutual development and benefits.



Product Description

Weiyou® ear tag applicator is made of high quality polyurethane (TPU), which is used for livestock and is for use with veterinary ear tags. Weiyou® ear tag applicator is a plier to install the ear tag on the pig, sheep, cattle etc. High quality, aluminum alloy material, shell high-grade painting material never rust, durable. Weiyou® ear tag applicator has resonable design according to the human plam, automatic lock , energy saving and easy to operate.

Specification

Product Name	Ear tag applicator
Property	animal ear tag plier
Usage	fixed the ear tag on the animal ear
Material	aluminum alloy
Size	24*6.5cm
Weight	340g
Consumable	ear tag pin/needle
Suitable	female ear tag, male ear tag
Package	OPP Bag, carton package for applicator

Features

A must for applying Global ear tags

Removable black plastic clip for applying electronic ear tags.



www.nbweiyou.com

Deep jaw makes proper placement easier.
Grip is designed to ease hand fatigue from repeated use.
each ear tag applicator with two super quality pin.
sandy surface for protecting slipping.
packaging: each piece with one blister packaging.
OEM packaging is available





www.nbweiyou.com



Ningbo Weiyou Import & Export Co., Ltd.
Phone:+86-13567833512
E-mail:_dario@nbweiyou.com



www.nbweiyou.com



Ningbo Weiyou Import & Export Co., Ltd.
Phone: +86-13567833512
E-mail: dario@nbweiyou.com





















